

## **NARKOMED AV2+ Ventilator**

### **TECHNICAL SERVICE BULLETIN # 171**

<b>Product:</b>	Narkomed Anesthesia Systems with AV2+ Ventilators (Narkomed 2C, Narkomed 4, Narkomed GS, Narkomed Mobile, Narkomed Military)
<b>Status:</b>	Class 3- Information
<b>Date:</b>	July 23, 2003
<b>Re:</b>	New Norgren Repair Kit
<b>Reference Doc:</b>	ECN 03-0065
<b>Reason/Rationale:</b>	To correct the occurrence of intermittent internal leakage across the regulator valve seat. This is observed as a continuous increase of pressure, as indicated on the Inspiratory Flow gauge, during expiration.
<b>Solution/Procedure:</b>	Install new Norgren Repair Kit with stronger spring.
<b>Devices affected:</b>	All AV2+ ventilators using Norgren Regulators (DMI Part number 4114252)
<b>Timing/When:</b>	Upon next service related to Inspiratory Flow Regulator symptom described above under "Reason/Rationale".
<b>Additional Info:</b>	<p>This part <b>replaces</b> Norgren Repair kit part number 4115151. Part number 4115151 will no longer be available.</p> <p>Narkomed MRI, MRI-2 Anesthesia Systems with AV2+ ventilators are excluded from using this repair part due to no serviceability below the top level assembly (to prevent damage to shielding).</p>
<b>Cost:</b>	Repair kit can be obtained at the then current list price.
<b>Ordering Info/Parts:</b>	New Part number: <b>4118241- Kit- Repair- Norgren Reg</b> This can be order thru your Customer Care Area Representative.
<b>Distribution:</b>	DrägerService, All Authorized Service Organizations

If you have any questions, please contact Technical Support by phone at 1-800-543-5047 or by fax at 1-215-723-5935.